



Scout 3 Compact Travel Scooter

HCPCS: K0800



INSET A




INSET B



INSET C



INSET D



INSET E



INSET F

- Quick and easy disassembly
- Interchangeable color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy free wheel operation
- Available with 12AH or 20AH batteries
- 20AH batteries increase scooter's maximum range

PRODUCTS	Description
SFSCOUT3	Red/Blue, Scout 3-Wheel Scooter Red/Blue, 1/ea
SFSCOUT3-EXT	Red/Blue, Scout 3-Wheel 20AH Battery Red/Blue, 1/ea

SPECIFICATIONS	Description
HCPCS	K0800
Weight Capacity	300 lbs.
Top Speed	4.25 mph
Maximum Range	9 miles (12AH) 15 miles (20AH)
Turning Radius	45.5"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions	42.5" (L) x 20.5" (W)
Floor-to-Seat Height	17" - 19"
Seat Dimensions	16.5" (W) x 16" (D)
Controller	CTE 45A
Motor	24V x 270W x 4700 rpm
Batteries (pair)	12V x 12AH, 12V x 20AH (EXT Model Only)
Charger	2A Offboard
Brakes	Electromagnetic
Front Wheels	8" x 2" Flat-free
Rear Wheels	8" x 2" Flat-free
Base Weight	49 lbs.
Battery Weight (pair)	20 lbs.

SPECIFICATIONS	Description
Seat Weight	16 lbs.
Total Weight	85 lbs.
Heaviest Piece Weight (Front Section)	29 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller	24 Months
Warranty on Drive Train Components	24 Months
Warranty on Batteries	12 Months
Wheel Base Dimensions (L x W)	32.7" x 16.5"

ACCESSORIES	Description
SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea
SF8050	Sun Shade, 1/ea

PARTS	Description
BC31001	
LRM412103-11	
BC31009N	
C11-076-00100	
SC31034	
BC31006	
BCKEY	
SC31008	
SC31009	
BC35051	
S35003-P	
SC31002	
SC31054	
SC31001	
SC31010	
C08-069-01000	
BC31021	
SC31038	
SC31043	
SC31044-D	
SC31046	
SC31039	
SC31024	
SC31027	
SC31019	
SC31041	
SC31041-D	
SC31042	
SC31042-D	
SC31040	
BC35045	
SC31035	
SC31035-SC	
BC35040R	
BC35040L	
BC953014	
S31052	
SC31023	
SC31018	