

Cobra GT4 Heavy Duty Scooter



- New, innovative design for a distinctive, stylish appearance
- One-touch, infinite tiller angle adjustment allows the user to easily position the tiller into the ideal driving position
- Ergonomic delta control handles makes it easier for users with limited strength and dexterity to operate the scooter
- Large, air-filled, low profile tires and full suspension for smoothest outdoor ride quality
- "Blade" style alloy wheel rims
- Efficient, extra bright LED lighting (Tail Lights and Turn Signals)
- Choice of Full size, Sport Style, Captain's Seats in 20" or 22" widths
- Seats are height adjustable, swivel and feature a fold-down, reclining backrest, and auto-style slides allow forward or backward adjustment
- Padded armrests feature easy turn cylinders for angle adjustment
- Positioning belt included
- "High/Low" switch for instant speed control
- Anti-tip wheels for greater safety
- Easy, freewheel mode

PRODUCTS	Description
COBRAGT420CS	1/cs
COBRAGT422CS	1/cs

SPECIFICATIONS	Description
Weight Capacity	450 lbs.
Top Speed	10 mph
Maximum Range	35 miles
Turning Radius	82"
Climbing Angle	8°
Ground Clearance	5.5"
Dimensions	63" (L) x 28.8" (W)
Floor-to-Seat Height	27" - 30"
Seat Dimensions	20" (W) x 20" (D) 22" (W) x 20" (D)
Controller	Dynamic Rhino II 160A
Motor	24V x 950W x 5200 rpm
Batteries (pair)	12V x 75AH
Charger	5A Offboard
Brakes	Electromagnetic
Front Wheels	14" x 14" Pneumatic
Rear Wheels	14" x 14" Pneumatic
Base Weight	191 lbs.
Battery Weight (pair)	94 lbs.
Seat Weight	55 lbs.

SPECIFICATIONS

Description

Total Weight	340 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronics	24 Months Limited
Warranty on Battery	12 Month Limited

ACCESSORIES

Description

ACCKIT-XL	1/cs
-----------	------

PARTS

Description

S08-062-00165	1/ea
C10-062-00200	1/ea
C10-051-00507	1/ea
C10-062-00300	1/ea
S12-051-00100	1/ea
C18-001-01501	1/ea
C37-001-02300	1/ea
C12-051-00200	1/ea
S08-062-00300	1/ea
S08-062-00477	1/ea
S37-062-00177	1/pr
S35-062-00100	1/ea
C35-045-00300	1/ea
S35-062-00300	1/ea
C35-062-00101	1/ea
S02-062-00100	1/ea
C02-029-00200	1/ea
C25-052-00100	1/ea
S02-001-00601	1/ea
S03-029-00200	1/ea
S03-062-00100	1/ea
C03-062-00300	1/ea
C03-022-02108	1/ea
C03-022-02208	1/ea
S05-062-00300	1/ea
C06-051-00100	1/ea
C06-051-00201	1/ea
C26-060-00501	1/ea
C13-062-00100	1/ea
C15-062-00102	1/ea
C15-062-00502	1/ea
C01-062-00300	1/ea
C08-062-00501	1/ea
S01-062-00100	1/ea
S08-062-00277	1/ea
S08-062-00777	1/ea
S07-062-00100	
S07-062-0020	
C07-016-00600	1/ea
C07-062-00200	1/ea
C07-005-00900	1/ea
C25-039-00100	1/ea
S08-062-00677	1/ea
C08-062-01600	1/ea
C25-022-00500	1/ea
C31-029-00101	1/ea