

# Cobra GT4 4-Wheel Heavy Duty Scooter

## Order Form/Price List



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>		<b>P.O. #</b>	<b>Date:</b>
<b>Bill to:</b>		<b>Ship to:</b>	

### Scooter with Seating

- Captain's Seat**
- COBRAGT420CS..... 20" (W) x 20" (D) .....\$8,645.00
  - COBRAGT422CS..... 22" (W) x 20" (D) .....\$9,009.00

### Batteries & Charger

- 12V x 75AH Batteries (pair) .....Included
- 5A Battery Charger .....Included

### Accessories

- AB1000 Armrest Bag .....\$25.00
- AB1120 Backpack .....\$56.00
  
- AH1000 Crutch/Cane Holder .....\$53.00
- AH2000 Cylinder Tank Carrier .....\$70.00
- AB2100 Rear Basket .....\$70.00
  
- AZ4000 Scooter Cover .....\$72.00
- ST301-COVER 20" Seat Cover .....\$118.17
- ST306-COVER 22" Seat Cover .....\$118.17
  
- AT1000 Scooter Trailer .....\$204.00
  
- AZ0060 Cup Holder .....\$30.00

### Specifications

HCPCS	N/A
<b>Performance</b>	
Maximum Weight Capacity	450 lbs.
Maximum Speed	10 mph
Maximum Range	35 miles
Maximum Grade/Incline	7.5°
Turning Radius	82"
<b>Dimensions</b>	
Overall Length	63"
Overall Width	25"
Floor to Seat Height	27"-30"
Ground Clearance	7"
<b>Power</b>	
Motor	24V x 950W x 5200rpm
Controller	Dynamic Rhino II 160A
Brakes	Electromagnetic
Freewheel Mode	Yes
<b>Tires (Pneumatic)</b>	
Front	14" x 4"
Rear	14" x 4"
<b>Weights</b>	
Base Weight	191 lbs.
Battery Weight	94 lbs.
Seat Weight	55 lbs.
Total Weight	340 lbs.

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.*

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.