


Cruiser II Standard Hemi Wheelchair

HCPCS: K0002



- Adjustable height, flip back, detachable desk length armrests
- 6" wheel lock extensions included
- Anti-tippers standard
- Seat belt included
- Built in seat depth extension easily adjusts from 16" – 18"
- Embossed vinyl upholstery is durable, lightweight, and easy to clean
- Carbon steel frame with silver vein finish
- Composite, mag-style wheels are lightweight and maintenance free
- 8" front casters
- Push-to-lock wheel locks
- Available with swing-away footrests or elevating legrests
- Dual axle provides multiple seat to floor height positions
- Configured in hemi position out of box for a lower seat to floor height

PRODUCTS	Description
K216ADDA-ELR	1/ea
K216ADDA-SF	1/ea
K218ADDA-ELR	1/ea
K218ADDA-SF	1/ea
K220ADDA-ELR	1/ea
K220ADDA-SF	1/ea

SPECIFICATIONS	Description
Width Open	(16", 18", 20") 23" 25" 27"
Width Closed	12"
Height	(Lowest Axle Position) 33"
Seat-to-Floor Height	17.5" Standard, Adjustable to 19.5"
Back Height	16"
Seat Depth	16" Standard, Adjustable to 18"
Length with Riggings	42"
Product Weight without Riggings	(16", 18", 20") 37 lbs. 39 lbs. 41 lbs.
Carton Shipping Weight	(16", 18", 20") 46 lbs. 50 lbs. 51 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS860	Universal Forward Stabilizer, 1/pr
STDSSHK-K3	Super Hemi-Kit, 1/ea
STDELR-AL	Articulating Legrest For Manual Wheelchairs, 1/cs

ACCESSORIES	Description
STDS820	Telescoping IV Pole Attachment, 1/ea
STDS1040S	Universal Cup Holder, 1/ea
STDS6008-1	Oxygen Cylinder Carry Bag for Wheelchair, 1/ea
14300	Tension Adjustable Back Cushion, 1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks, 1/ea
STDS834	Overhead Anti-Theft Device Single Pole, 1/ea
STDS806	Anti-Fold Bar, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea
STDS5050	Easy-To-Clean Surface. Accommodates Wheelchairs, 1/ea
WASR	Limb Support, Right, 1/ea
WASL	Limb Support, Left, 1/ea
STDS835	Chart Carry Pocket, 1/ea
STDS836	Chart Carry Pocket, 1/ea