

Titan AXS Mid-Wheel Drive Power Wheelchair Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Powerchair with Seating

High Back Captain's Seat

- TITANAXS-1616CS16" (W) x 16" (D)\$6,725.50
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.
- TITANAXS-1618CS16" (W) x 18" (D)\$6,725.50
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.
- TITANAXS-1816CS18" (W) x 16" (D)\$6,725.50
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.
- TITANAXS-18CS18" (W) x 18" (D)\$6,725.50
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.
- TITANAXS-20CS20" (W) x 18" (D)\$6,891.50
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.
- TITANAXS-2020CS20" (W) x 20" (D)\$6,957.00
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.
- TITANAXS-22CS22" (W) x 20" (D)\$7,056.00
HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.

Standard Armrests

- ARM-L.....Included
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests
- ARM-RIncluded
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests

Controller Arm

- Adjustable Length.....Included
- AA4100 – Swing-Away Controller Arm (HCPCS E1028)\$188.37

Front Riggings

- Adjustable FootplateIncluded
- AE2500/LK3JELR – Elevating Legrest Kit\$261.89
HCPCS: E0990, Includes swing-away elevating legrest with attachment bracket

Batteries & Charger

- U1BATTERY (pair)Included
HCPCS: E2365, 12V x 35AH
- 4A Battery ChargerIncluded

Accessories

- STDS850-PW Positioning/Safety Belt (HCPCS E0978).....Included
- SF8030 Crutch/Cane Holder (HCPCS E2207)\$101.23
- SF8010 Cylinder Tank Carrier (HCPCS E2208)\$150.91

Specifications

HCPCS

TITANAXS-1616CS	K0823
TITANAXS-1618CS	K0823
TITANAXS-1816CS	K0823
TITANAXS-18CS	K0823
TITANAXS-20CS	K0823
TITANAXS-2020CS	K0823
TITANAXS-22CS	K0823

Performance

Maximum Weight Capacity	300 lbs.
Maximum Speed	4 mph
Maximum Range	19 miles
Maximum Grade/Incline	6°
Turning Radius	20"

Dimensions

Overall Length	39.25"
Overall Width	23.5"
Floor to Seat Height	20.75 – 21.75"
Ground Clearance	2.5"

Power

Motor	Two 24V x 320W x 4600rpm
Controller	Dynamic Linx 40A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires

Casters	6" x 2" Solid
Drive Wheels	10" x 3" Flat Free

Weights

Base Weight	87 lbs.
Battery Weight	44 lbs.
Seat Weight	40 lbs.
Total Weight	171 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.