


Titan AXS Mid-Wheel Drive Powerchair



- Captain's seat, available in 7 different width and depth combinations, is height adjustable, reclines, swivels in 45° increments, and has a height adjustable headrest
- Padded, contoured armrests with height, width and angle adjustment
- Adjustable length controller mount
- Interchangeable color panels in red and blue included
- Flip-up footplate with height, depth and angle adjustment
- Positioning belt included
- Programmable controller with Bluetooth electronics for easy trouble-shooting and service
- Flat-free, non-marking tires
- Easy freewheel operation

PRODUCTS

Description

TITANAXS-1616CS	K0823, 1/ea
TITANAXS-1618CS	K0823, 1/ea
TITANAXS-1816CS	K0823, 1/ea
TITANAXS-18CS	K0823, 1/ea
TITANAXS-20CS	K0823, 1/ea
TITANAXS-2020CS	K0823, 1/ea
TITANAXS-22CS	K0823, 1/ea

SPECIFICATIONS

Description

Weight Capacity	300 lbs.
Top Speed	4 mph
Estimated Range	19 miles
Max. Climbing Angle	6°
Turning Radius	20"
Ground Clearance	2.5"
Overall Dimensions	39.25" (L) x 23.5" (W) x 44" (H)
Seat Dimensions	16" 18" 20" 22"
Depth	16" 18" 20"
Back Height	18"
Seat-to-Floor Height	20.75" – 21.75"
Motor	320W x 4600RPM
Controller	Dynamic Linx 40A
Batteries	12V x 33AH
Battery Charger	4A
Brakes	Electromagnetic
Freewheel Mode	Yes

SPECIFICATIONS	Description
Tires	Flat Free
Caster Wheels	6" x 2"
Drive Wheels	10" x 3"
Base Weight	87 lbs.
Battery Weight (Pair)	44 lbs.
Seat Weight	40 lbs.
Total Weight	171 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronics	24 Months
Warranty on Batteries	12 Months

ACCESSORIES	Description
AE2500	ELR Bracket, 1/ea
LK3JELR	ELR Set, 1/pr
AA4100	Swing-Away Joystick, 1/ea
SF8010	Oxygen Tank Holder, 1/ea
SF8030	Cane/Crutch Holder, 1/ea
ST205-COVER	18" Seat Cover, 1/ea
ST301-COVER	20" Seat Cover, 1/ea
ST306-COVER	22" Seat Cover, 1/ea