



Phoenix HD 3 Heavy-Duty Travel Scooter

HCPCS: K0801



INSET A



INSET B



INSET C



INSET D



INSET E

- Quick and easy disassembly
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with two-tone upholstery
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Swivel seat with fold-down backrest
- Easy free wheel operation

PRODUCTS

Description

PHOENIXHD3	Phoenix HD 3Whl Travel POV 17.5 Seat, 1/cs
PHOENIXHD3-20	Phoenix HD 3Whl Travel POV 20 Seat, 1/cs

SPECIFICATIONS

Description

HCPCS	K0801
Weight Capacity	350 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	43"
Climbing Angle	6°
Ground Clearance	3.75"
Dimensions	39" (L) x 22.5" (W)
Floor-to-Seat Height	20" - 22"
Drive System	24(V) x 350(W) x 3600 rpm
Seat Type	Folding
Seat Dimensions	17.5" (W) x 17" (D)
Motor	24V x 350W x 3600 rpm
Batteries (pair)	12V x 20AH
Charger	1.5A Offboard
Brakes	Electromagnetic
Front Wheel	9"x 3" Flat-free
Rear Wheels	9"x 3" Flat-free
Base Weight	63 lbs.

SPECIFICATIONS	Description
Battery Weight (pair)	31 lbs.
Seat Weight	21 lbs.
Total Weight	115 lbs.
Heaviest Piece Weight	32 lbs. (Front Section)
Warranty on Frame	Lifetime
Warranty on Electronic Controller	24 Months
Warranty on Drive Train Components	24 Months
Wheel Base Dimensions (L x W)	29.1" x 19.1"
Warranty on Batteries	12 Months

ACCESSORIES	Description
SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea
SF8050	Sun Shade, 1/ea

PARTS	Description
S31009	1/ea
S31003	1/ea
C20-007-0020	
S31004GT	
S31209BK	1/ea
S31006	1/ea
S3KEY	
S31022	1/ea
S31001	1/ea
LRM412103-11	1/ea
S31012	1/ea
S31122	1/ea
S31208B	1/ea
S31147	1/ea
BC35051	1/ea
S350109	1/ea
S31028	1/ea
S31021	1/ea
S31036-3HD	1/ea
S350127	1/ea
S31012B	
S31012R	
S31142-3HD	1/ea
9503H79021	1/pr
S31141-3	
S31011	1/ea
S350169715	1/ea
S350169716	1/ea
S350169717	1/ea
S350110	1/ea
S350113	1/ea
S350114	1/ea
C08-035-01500	
S350169707	1/ea
S350169706-HD	1/ea
S350169706 -HDCP	
LRM412123	1/ea
LRM302111	1/ea