



**Ventura 4  
Midsize Scooter**

HCPCS: K0807



INSET A



INSET B



INSET C



INSET D



- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Seat position can be adjusted forward or backward
- Front and rear suspension
- Programmable controller
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Delta tiller
- Easy free wheel operation

PRODUCTS	Description
VENTURA418FS	VENTURA 4-WHEEL 18 FOLDING SEAT RED/BL, 1/cs

SPECIFICATIONS	Description
HCPCS	K0807
Weight Capacity	400 lbs.
Top Speed	6 mph
Maximum Range	16 miles
Turning Radius	54"
Climbing Angle	7.5°
Ground Clearance	2.5"
Dimensions	48" (L) x 25" (W)
Floor-to-Seat Height	22" - 24"
Seat Type	Folding
Seat Dimensions	18" (W) x 17.5" (D)
Motor	24V x 350W x 5100 rpm
Batteries (pair)	12V x 33AH
Charger	4A Offboard
Brakes	Electromagnetic
Rear Wheels	10" x 3" Flat-free
Base Weight	99 lbs.
Battery Weight (pair)	28 lbs.

SPECIFICATIONS	Description
Heaviest Piece Weight	70 lbs.
Total Weight	153 lbs.
Seat Weight	26 lbs.
Front Wheels	10" x 3" Flat-free
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Batteries	12 Months

ACCESSORIES	Description
1/ea	SF8010
1/ea	SF8020
1/ea	SF8030
1/ea	SF8050

PARTS	Description
VENT-01	
VENT-04	
VENT-05	
S31001	
LRM412103-11	
S31208B	
VENT-72	
S31122	
GT8072W38	
GT8072W43	
GT8072W41	
GT8072W42	
S31012	
GT8072W32	
GT8072W44	
VENT-76	
S41160	
VENT-03	
VENT-53	
VENT-54	
GT8072W35	
GT8072W36	
GT8072W53	
S35003-PL	
VENT-55	
VENT-61	
VENT-62	
VENT-06	
VENT-02	
VENT-07	
S31147	
VENT-44	
VENT-10	
VENT-27	
VENT-23	
VENT-26	
VENT-24	
VENT-25	
VENT-56	
VENT-59	