


Ventura 3 DLX Midsize Scooter



INSET A



INSET B



INSET C



INSET D



INSET E



INSET F

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- 33 AH batteries for extended range
- Floorboard with 23" of foot space
- Ergonomic throttle control
- Height-adjustable swivel seat features fold-down backrest and adjustable padded armrests
- Seat position can be adjusted forward or backward
- Rear suspension and Anti-tip wheels
- Programmable controller
- Seat back reclines and features adjustable headrest
- Available with 18" or 20" width Captain's Seat
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, Turn signals/hazard lights and Rearview mirrors
- Easy free wheel operation

PRODUCTS

PRODUCTS	Description
VENTURA318CS	VENTURA 3-WHL DELUXE 18CAPT.SEAT RED/BL, 1/cs
VENTURA320CS	VENTURA 3-WHL DELUXE 20CAPT.SEAT RED/BL, 1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
HCPCS	N/A
Weight Capacity	400 lbs.
Top Speed	6 mph
Maximum Range	16 miles
Turning Radius	53"
Climbing Angle	7.5°
Ground Clearance	2.5"
Dimensions	48" (L) x 25" (W)
Floor-to-Seat Height	20.5" - 23.5"
Seat Type	Captain's
Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)
Motor	24V x 350W x 5100 rpm
Batteries (pair)	12V x 33AH
Charger	4A Offboard
Brakes	Electromagnetic
Front Wheels	10" x 3" Flat-free

SPECIFICATIONS	Description
Rear Wheels	10" x 3" Flat-free
Base Weight	88 lbs.
Battery Weight (pair)	46 lbs.
Seat Weight	40 lbs.
Heaviest Piece Weight	60 lbs.
Total Weight	174 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Batteries	12 Months

ACCESSORIES	Description
SF8010	1/ea
SF8020	1/ea
SF8030	1/ea
SF8050	1/ea

PARTS	Description
VENT-04	
VENT-05	
VENT-01	
S31208B	
S31122	
GT8072W38	
GT8072W43	
S31012	
S31001	
LRM412103-11	
GT8072W32	
GT8072W44	
GT8072W41	
GT8072W42	
VENT-72	
S341160	
VENT-54	
VENT-53	
VENT-03	
VENT-61	
VENT-62	
VENT-06	
VENT-07	
GT8072W35	
GT8072W36	
GT8072W53	
VENT-45	
VENT-46	
VENT-02	
VENT-13	
VENT-48	
VENT-49	
VENT-58	
VENT-59	
VENT-60	
VENT-50	
VENT-44	
S35003-PL	