

Pediatric Kanga TS Inclusive System

HCPCS: E1236



FIGURE A



- Durable lightweight aluminum frame
- Built-in growth frame is width adjustable and can be configured for 10", 12" and 14" Seating system and positioning accessories must be changed accordingly
- Height- and angle-adjustable flip-back armrests
- Height-adjustable push handles
- Tilt-in-Space up to 45 degrees
- 6" of adjustable seat depth
- Seat-to-back angle adjustable
- Angled seat to prevent forward sliding
- Swing-away, height-adjustable footrests with foot and ankle positioner
- Angle-adjustable footplates
- 20" rear semi-pneumatic quick-release wheels
- Frame accommodates most custom rehab seating and positioning systems
- Curb-assist levers
- Height-adjustable anti-tippers
- Crash tested and meets bus transit requirements
- Folds for easy transport (Figure A)

PRODUCTS

Description

KG 1000	Red, 1/bx
KG 1400	Red, 1/bx

SPECIFICATIONS

Description

Push Handle Height	38.75" - 43.5"
Armrest Height	10" - 12"
Back Cane Height	18"
Back Angle Adjustment	5° - 25°
Carton Shipping Weight	53 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

KG 8024	Kanga TS Upper Extremity Tray, 1/bx
KG 8026	Kanga TS Canopy, 1/bx
KG 8000	Headrest for Kanga TS Folding Wheelchairs, 1/bx
KG 8045	Heel and Toe Loops for Pediatric Kanga TS Wheelchair, 1/pr
KG 1040	Adductors for Kanga TS Folding Wheelchairs, 1/pr
KG 1041	Abductors for Kanga TS Folding Wheelchairs, 1/pr
KG PELR	Elevating Leg Rests, 1/pr
KG AELR	Elevating Leg Rests, 1/pr